

# SHENYUN

## STERILE SINGLE USE NEGATIVE PRESSURE



## ■ GUIDE NEEDLE



- Made of stainless steel.
- Sharp, smooth and burr free.
- Protection cover.
- Customized specification is available upon request.

Code(Not foldable)	Size	Needle diameter/mm	length/mm
SY03-08-01	10Fr	3.2	128
SY03-08-02	12Fr	4.0	128
SY03-08-03	14Fr	4.8	128
SY03-08-04	16Fr	5.3	128
SY03-08-05	18Fr	6.0	128
SY03-08-06	19Fr	6.3	128
SY03-08-07	20Fr	6.7	128
SY03-08-08	22Fr	7.3	128
SY03-08-09	24Fr	8.0	128
Code (foldable)	Size	Needle diameter/mm	length/mm
SY03-08-01	10Fr	3.2	145
SY03-08-02	12Fr	4.0	145
SY03-08-03	14Fr	4.8	145
SY03-08-04	16Fr	5.3	145
SY03-08-05	18Fr	6.0	145
SY03-08-06	19Fr	6.3	145
SY03-08-07	20Fr	6.7	145
SY03-08-08	22Fr	7.3	145
SY03-08-09	24Fr	8.0	145



# PRODUCT TRADING

## STERILE SINGLE USE NEGATIVE PRESSURE

The disposable surgical drainage set is a common medical device used for drainage of body fluids after surgery or in trauma treatment.

Its design features and components are designed to provide an effective, sterile drainage solution while being easy to use and replace.

### ■ Features

- **Disposable:** Sterile before use, discarded after use to prevent cross infection.
- **Efficient drainage:** The porous tube improves drainage efficiency and ensures smooth discharge of body fluids.
- **Safety:** Biocompatible materials: Avoid irritation and allergic reactions to tissues.
- **Comfort:** Lightweight design makes it easy for patients to carry and move around.

### ■ Precautions

- **Aseptic operation:** Maintain aseptic operation during use to avoid contamination of the drainage tube and drainage bag.
- **Firm fixation:** Ensure that the drainage tube is firmly fixed to avoid displacement or falling off.
- **Replace in time:** Replace the drainage bag regularly to prevent excessive fluid accumulation leading to reflux or infection.